




TYPE III DR-92 G3500

GM 3500 6.6 DURAMAX TURBO DIESEL

CHASSIS SPECIFICATIONS

Chassis	GM 3500 6.6 Duramax 6600 Turbo Diesel
Chassis, Alternator	OEM Package
Chassis, Color	Interior Gray
Chassis, Floor Covering	OEM
Chassis, Mirrors	Remote Control with Convex
Chassis, Seats	Cloth Captains Chairs
Chassis, Throttle	OEM system
Chassis, Trim	High Trim Level



CONVERSION SPECIFICATIONS

Conversion	TraumaHawk, High Trim
Module	145" x 92" x 66" Interior Headroom
Attendant Seat	Captain's Chair on Storage Box
Batteries	(1) Underhood and (1) in RF Compartment
Battery Switch	Cole Hersee 2484-16
Backboard Storage	Left Rear Compartment
Back Up Alarm	With Auto Reset Switch
Biowaste System	Curbside, at Head of Squad Bench
Bumper, Rear	Treadplate, with Center Flip Up Section
Cabinet, Action Area	Mica Covered Wood
Cabinet, Behind Attendant Seat	With Outside Access for Spare Tire
Cabinet, Electrical	Above Walkthru with Lever Latch
Cabinet, Laminate	Choice of Optional Colors
Cabinet, Partition Door	Sliding with Sliding Window
Cabinet, Polycarbonate	Gray Tinted or Clear with Full Handles
Cabinet, RF	Open for Easy Inside Access
Cabinet, Storage	Door Under Action Area for Electrical Components
Cabinet, Straps	(1) in Left Rear Cabinet for Stairchair
Cabinet, Streetside	Large Openings with Set Back Dividers
Cabinet, Tambour Doors	(1) Total in Patient Area
Coax Cables	(1) Tagged at Both Ends
Compartment, Lining	Diamond Plate
Console	On Engine Cover
Cot Mount	Stryker
Door Hold Opens	Grabbers on Rear Doors
DOT Package	Shipped Loose
Drip Rails	Bright Finish Aluminum
Electrical System	With Custom Electrical Schematics
Electrical Wiring	Color, Name & Number Coded
Fasteners, Exterior	AVK for Attachment of Door Panels
Fender Flares	Rubber
Flooring	Choice of Optional Colors
Fuel Fill	Aluminum Housing
Ground Straps	(4) with Star Washers
Handles, Interior	Handicap Style on all Entry Doors

Headliner	(1) Piece Front to Rear, Gloss White
Heat-AC	Ceiling Ducted
Insulation	Reflectix Circumferential Package
IV Hooks	(2) Perko Hooks Installed Overhead
KKK Package	To Present Revision
License Plate Holder	Recessed in Rear Kickplate
Lights, Action Area	12vdc LED
Lights, Clearance	LED On Front & Rear
Lights, Dome	(6) Total LED, (3) Curbside (3) Streetside
Lights, Indicator	Door, Compartment Ajar, & Battery on Cab Console
Lights, Load	(2) Whelen on Rear
Lights, Scene	(2) Whelen on Each Side
Lights, Side Marker	On Sides of Body
Lights, Tail	Triple Light Cluster on Rear
Lights, Third Brake	Mounted over Rear Doors
Lights, Warning Package	Whelen Perimeter Lighting
Outlets, 125vac	(2) Total, (1) Action Area (1) RF-ALS
Outlets, 12vdc	(2) Cigar Outlets
Overhead Grabrail	72" with Rounded Ends, Recessed
Oxygen Cylinder Rack	Aluminum Rack, LF Compartment
Oxygen Outlets	(2) Total, Ohio Style
Paint	OEM White Industry Match
Exterior Striping	Single Belt
Patient Floor	Plywood, Full Length & Width, Insulated
Reflectors	(6) Total, (4) Rear & (2) Front
Shoreline	15 Amp Inlet on Streetside
Siren	Whelen WS-295
Siren Speaker	Dual 100watt
Skirtrails	Diamond Plate
Spare Tire	Mounted in Left Front Compartment
Suction System	SSCOR with Disposable Container, Electric Pump
Undercoat	Per Industry Guidelines
Upholstery	Choice of Optional Colors
Windows, C/S Upper	Std. Tint Fixed
Windows, Rear Upper	Std. Tint Fixed



*Base conversion shown—additional configurations and options available as part of AEV's ongoing program of product improvement. Specifications are subject to change without notice.

165 AMERICAN WAY, JEFFERSON, NC 28640 | 336.982.9824 | 800.374.9749 | FAX: 336.982.9826 | WWW.AEV.COM

TYPE III DR-92 G3500 6.6 DURAMAX TURBO DIESEL V10.1





GM 3500 6.6 DURAMAX TURBO DIESEL

TYPE III DR-92 G3500



TYPE III DR-92 G3500

GM 3500 6.6 DURAMAX TURBO DIESEL