



TYPE III DR-92 E-350

CHASSIS SPECIFICATIONS

Chassis.....	Ford E-350
Chassis, Alternator.....	OEM Package
Chassis, Ambulance.....	47A Ford Ambulance Package
Chassis, Brakes.....	4-Wheel Anti Lock
Chassis, Color Interior.....	Gray
Chassis, Floor Covering.....	OEM, rubber mat
Chassis, Headlights.....	Euro Style with Daytime Running Lights
Chassis, Mirrors.....	OEM Power with OEM Convex
Chassis, Power.....	Windows and Door Locks
Chassis, Radio.....	AM/FM/CD, Stereo with Aux Input
Chassis, Seats.....	Cloth Captains Chairs with Inboard Armrests
Chassis, Throttle.....	Activated thru parking brake
Chassis, Transmission.....	Automatic with Overdrive
Chassis, Trim.....	XL, High Trim Level
Chassis, Wheelbase.....	138"



CONVERSION SPECIFICATIONS

Conversion.....	TraumaHawk, High Trim
Module.....	145" x 92" x 66" Interior Headroom

Attendant Seat.....	Captains Chair on Storage Box
Batteries.....	(2) Total, OEM Locations
Battery Switch.....	Cole Hersee 2484-16
Backboard Storage.....	Left Rear Compartment
Back Up Alarm.....	With Auto Reset Switch
Biowaste System.....	Curbside, at Head of Squad Bench
Bumper, Rear.....	Treadplate, with Center Flip Up Section
Cabinet, Action Area.....	Mica Covered Wood
Cabinet, Behind Attendant Seat..	With Outside Access for Spare Tire or Stairchair
Cabinet, Electrical.....	Above Walkthru with Lever Latch
Cabinet, Laminate.....	Choice of Optional Colors
Cabinet, Partition Door.....	Sliding with Sliding Window
Cabinet, Polycarbonate.....	Gray Tinted or Clear with Full Handles
Cabinet, RF.....	Open for Easy Inside Access
Cabinet, Storage.....	Access Under Action Area for Electrical Components
Cabinet, Streetside.....	Large Openings with Set Back Dividers
Cabinet, Tambour Doors.....	(1) Total in Patient Area
Coax Cables.....	(1) Tagged at Both Ends
Compartment, Lining.....	Diamond Plate
Console.....	ABS on Ford Engine Cover
Cot Mount.....	Stryker
Door Hold Opens.....	Grabbers on Rear Doors
DOT Package.....	Shipped Loose
Drip Rails.....	Bright Finish Aluminum
Electrical System.....	With Custom Electrical Schematics
Electrical Wiring.....	Color, Name & Number Coded
Fasteners, Exterior.....	AVK for Attachment of Door Panels
Fender Flares.....	Rubber
Flooring.....	Choice of Optional Colors
Fuel Fill.....	Aluminum Housing
Ground Straps.....	(4) with Star Washers
Handles, Interior.....	Handicap Style on all Entry Doors

Headliner.....	(1) Piece Front to Rear, Gloss White
Heat-AC.....	Ceiling Ducted
Insulation.....	Reflectix Circumferential Package
IV Hooks.....	(2) Perko Hooks Installed Overhead
License Plate Holder.....	Recessed in Rear Kickplate
Lights, Action Area.....	12vdc LED with Switch
Lights, Clearance.....	LED On Front & Rear
Lights, Dome.....	(6) Total LED, (3) Curbside (3) Streetside
Lights, Indicator.....	Door, Compartment Ajar, & Battery on Cab Console
Lights, Load.....	(2) Whelen on Rear
Lights, Mica.....	(2) Whelen on Each Side
Lights, Side Marker.....	On Sides of Body
Lights, Tail.....	Triple Light Cluster on Rear
Lights, Third Brake.....	LED Mounted over Rear Doors
Lights, Warning Package.....	Whelen Perimeter Lighting
Outlets, 125vac.....	(2) Total, (1) Action Area and (1) RF-ALS
Outlets, 12vdc.....	(2) Cigar Outlets
Overhead Grabrail.....	72" with Rounded Ends, Recessed
Oxygen Cylinder Rack.....	Aluminum Rack, LF Compartment
Oxygen Outlets.....	(2) Total, Ohio Style
Paint.....	OEM White Industry Match
Exterior Striping.....	Single Belt
Patient Floor.....	Plywood, Full Length & Width, Insulated
Reflectors.....	(6) Total, (4) Rear & (2) Front
Shoreline.....	15 Amp Inlet on Streetside
Siren.....	Whelen WS-295
Siren Speaker.....	Dual 100watt
Skirtrails.....	Diamond Plate
Spare Tire.....	Mounted in Left Front Compartment
Straps.....	(2) Total, (1) Stairchair and (1) Backboards
Suction System.....	SSCOR with Disposable Container, Electric Pump
Undercoat.....	Per Ford QVM Guidelines
Upholstery.....	Choice of Optional Colors
Windows, C/S Upper.....	Std. Tint Fixed
Windows, Rear Upper.....	Std. Tint Fixed



*Base conversion listed—additional configurations and options available as part of AEV's ongoing program of product improvement. Specifications are subject to change without notice.

165 AMERICAN WAY, JEFFERSON, NC 28640 | 336.982.9824 | 800.374.9749 | FAX: 336.982.9826 | WWW.AEV.COM

FORD DR-92 E-350 V11.1





TYPE III DR-92 E-350



TYPE III DR-92 E-350