



# TYPE I D4500

DODGE 4500/5500, 6.7L TURBO DIESEL



**SAFETY IS NOT AN OPTION**



AMERICAN EMERGENCY VEHICLES



# TYPE I D4500

DODGE 4500/5500, 6.7L TURBO DIESEL

**AEVE**  
AMERICAN EMERGENCY VEHICLES



# TYPE I D4500

DODGE 4500/5500, 6.7L TURBO DIESEL

## CHASSIS SPECIFICATIONS

Chassis . . . . .	Dodge 4500/5500, 6.7L Turbo Diesel
Chassis, Alternator . . . . .	OEM
Chassis, Color Interior . . . . .	Gray
Chassis, Floor Covering . . . . .	OEM
Chassis, Mirrors . . . . .	OEM
Chassis, Seats . . . . .	Cloth Captains Chairs
Chassis, Throttle . . . . .	OEM system
Chassis, Trim . . . . .	SLE, High Trim Level



## CONVERSION SPECIFICATIONS

Conversion . . . . .	TraumaHawk, High Trim
Module . . . . .	148" x 95" x 68" Interior Headroom

Attendant Seat . . . . .	Captains Chair on Storage Box
Batteries . . . . .	(2) Underhood
Battery Switch . . . . .	Cole Hersee 2484-16
Backboard Storage . . . . .	Right Rear Compartment, Vertical
Back Up Alarm . . . . .	With Auto Reset Switch
Biowaste System . . . . .	Streetside
Bumper, Rear . . . . .	Treadplate, with Center Flip Up Section
Cabinet, Action Area . . . . .	Laminated Countertop
Cabinet, CPR Seat . . . . .	Streetside with Telemetry Area
Cabinet, Electrical . . . . .	Behind Attendants Seat
Cabinet, Laminate . . . . .	Choice of Optional Colors
Cabinet, Polycarbonate . . . . .	Gray Tinted or Clear with Full Handles
Cabinet, RF . . . . .	Wood Doors with Inside/Outside Access
Cabinet, Tambour Doors . . . . .	(1) Total in Patient Area
Coax Cables . . . . .	(1) Tagged at Both Ends
Compartment, Lining . . . . .	Diamond Plate
Compartment, Shelves . . . . .	(1) Total Aluminum in Exterior Compartment
Compartment, Straps . . . . .	(1) Total in Right Rear Compartment for Backboards
Console . . . . .	Wood Between Cab Seats
Cot Mount . . . . .	Stryker
Door Hold Opens . . . . .	Grabbers on Rear Doors
DOT Package . . . . .	Shipped Loose
Drip Rails . . . . .	Bright Finish Aluminum
Electrical System . . . . .	With Custom Electrical Schematics
Electrical Wiring . . . . .	Color, Name & Number Coded
Fasteners, Exterior . . . . .	AVK for Attachment of Door Panels
Fender Flares . . . . .	Aluminum
Flooring . . . . .	Choice of Optional Colors
Flooring, Roll Up . . . . .	3" Recessed
Fuel Fill . . . . .	Aluminum Housing
Ground Straps . . . . .	(4) with Star Washers

Handles, Interior . . . . .	Handicap Style on all Entry Doors
Headliner . . . . .	(1) Piece Front to Rear, Gloss White
Heat-AC . . . . .	Ceiling Ducted
Insulation . . . . .	Reflectix Circumferential Package
IV Hooks . . . . .	(2) Perko Hooks Installed Overhead
KKK Package . . . . .	To Present Revision
License Plate Holder . . . . .	Recessed in Rear Kickplate
Lights, Action Area . . . . .	12vdc LED
Lights, Clearance . . . . .	LED On Front & Rear
Lights, Dome . . . . .	(6) Total, (3) Streetside and (3) Curbside
Lights, Indicator . . . . .	Door, Compartment Ajar, & Battery on Cab Console
Lights, Load . . . . .	(2) Whelen on Rear
Lights, Scene . . . . .	(2) Whelen on Each Side
Lights, Side Marker . . . . .	On Sides of Body
Lights, Tail . . . . .	Triple Light Cluster on Rear
Lights, Third Brake . . . . .	Mounted over Rear Doors
Lights, Warning Package . . . . .	Whelen Perimeter Lighting
Outlets, 125vac . . . . .	(2) Total, (1) A/A and (1) RF-ALS
Outlets, 12vdc . . . . .	(2) Cigar Outlets
Overhead Grabrail . . . . .	72" with Rounded Ends, Recessed
Oxygen Cylinder Rack . . . . .	Aluminum Rack, LF Compartment
Oxygen Outlets . . . . .	(2) Total, Ohio Style
Paint . . . . .	OEM White Industry Match
Exterior Striping . . . . .	Single Belt
Patient Floor . . . . .	Plywood, Full Length & Width, Insulated
Reflectors . . . . .	(6) Total, (4) Rear and (2) Front
Shoreline . . . . .	15 Amp Inlet on Streetside
Siren . . . . .	Whelen WS-295-HFRS
Siren Speaker . . . . .	Dual 100 watt
Skirtrails . . . . .	Diamond Plate
Spare Tire . . . . .	Shipped Loose
Suction System . . . . .	SSCOR with Disposable Container, Electric Pump
Undercoat . . . . .	Per QVM Guidelines
Upholstery . . . . .	Choice of Optional Colors
Windows, C/S Upper . . . . .	Std. Tint Fixed
Windows, Pass Thru . . . . .	Sliding Lexan with Bellows
Windows, Rear Upper . . . . .	Std. Tint Fixed



\*Base conversion shown—additional configurations and options available as part of AEV's ongoing program of product improvement. Specifications are subject to change without notice.