

Introducing the Safest Ambulance on the Streets

AMR's most advanced "next generation C-3 ambulance" makes medical transportation safer for medics, in turn allowing for better patient care both before & during transport.

- High output Heat/AC system accommodates variable filtration capabilities, includes HVAC ducted through patient compartment headliner. HEPA filtration continuously cleans the air.
- Full Heat and A/C controls from either patient area switch panel and cab console.

- P6 6-point seat belt systems on both curbside seating locations with built in seat sensors for driving identification of seat belted locations.

- Special lighting package includes amber LED caution lights that flash at motorist eye level, as well as low amperage LED emergency warning lights

- 3 point integral seat belt attendant's seat with captain's chair including an internal fold down child safety seat for transport convenience of small children.



- Easy to identify overhead grab rail and restraint netting in patient area.

- 3rd battery for emergency start backup system for added security.

- Dual full fuction onboard suction systems for use from multiple seating locations.

- Brake and Turn indicator lights in patient area over rear doors to warn attendant of pending turn or stop.

- Low amperage LED side scene and rear loading lights to help maintain chassis electrical system integrity.

- Fluorescent lights in patient area center pad with timer activated by the opening of the curbside entry door to allow preset activation with auto shutoff programmed up to 1 hour for vehicle checkout with rear master switch turned off.

- Front body LED turn arrows to warn oncoming traffic of pending turn in emergency situation.



- Fly by wire, multiplexed electrical system with State of the Art technology includes the capability for remote programming updates from the factory to update programming, trouble shoot for change switch activations and features.

- Cabinet at the foot of the squad bench for crew required equipment with easy access without reaching over the patient during transport.



- New low amperage LED patient area dome lights with added intensity control and longer life expectancy.

- On board camera system with up to four cameras for view of patient area, backup, right fender lane change and forward for recording capability.

- Micro 4 backup sensing device for assistance with backing the vehicle with audible cue to the driver with number of feet to obstruction.

- High conspicuity graphics package for both day and night visibility on the road and at intersections.

- Programmable signs on streetside and curbside of body for notification during both transport and emergency transports.

- Seating includes vacuum formed upholstery for a seamless joint to allow for easy and constant cleaning like in a hospital environment.

- Crash stable LifePak bracket located in the action area with spring loaded release for easy access for use on scene and during transport.



New features added to Concept Ambulance “C-3” based on feedback from EMS professionals across the nation:



- Digital O2 Contents Monitor
- Seat sensing system for squad bench
- Up to four camera system with recording
- HEPA Air Filtration
- High Conspicuity Graphics
- Multiplex Programmable Electrical System
- Multiple Patient Area Switch Panels
- P6 Seating for squad bench
- Load Manager/Sequential Lighting
- Crash Stable LifePak bracket

- Safety Net at Head of squad bench
- Easy clean interior
- Ducted Heat and A/C System
- Low amperage electrical components
- Integrated Child Safety Seat
- Variable Speed Heat and A/C fan control
- Ergonomic Switch Panel in Cab
- Bio Waste system in Action Area

